




Test Certificate

Part 1 Liter Vessel  
Part Number 14143-01  
Serial Number 104301 - 104515



		USP / Ph.Eur. Tolerances	
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm
		Height b	200 mm
		With hemispherical bottom	160 – 210 mm
Material	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2025-11-11

Buner Cornelius  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland