




Part 1 Liter Vessel
 Part Number 14143-11
 Serial Number P35 000800 – P35 000907



		USP / Ph.Eur. Tolerances	
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm
		Height b	200 mm
		With hemispherical bottom	160 – 210 mm
Material	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 26.03.2021

26-03-2021 11:11:02

Moreshwar Pulekar
 SOTAX Quality Management
 EL 142/143 | Ground Floor | TTC Industrial area | MIDC Mahape | Navi Mumbai 400701 | India

Test Certificate