




Part 1 Liter Vessel
 Part Number 14143-11
 Serial Number P35 001529 – P35 001684



		USP / Ph.Eur. Tolerances	
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm
		Height b	168.2 mm
		With hemispherical bottom	160 – 210 mm
Material	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2022-06-27

27-06-2022 14:09:21

Moreshwar Pulekar
 SOTAX Quality Management
 C-64 | First Floor | TTC Industrial area | Turbhe | Navi Mumbai 400705 | India

Test Certificate